

Product datasheet for **HP214986**

IQCA1 Human qPCR Primer Pair (NM_024726)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79781
Forward Sequence:	CTACTCTTCGCCAGGTGTCCAT
Reverse Sequence:	GGAGCTTTCTCATGCACTGCTG
ACCN:	<u>BC028699</u> , <u>NM_024726</u> , <u>NM_024726.1</u> , <u>NM_024726.2</u> , <u>NM_024726.3</u> , <u>NM_024726.4</u> , <u>BC028699.1</u> , <u>BC043504</u> , <u>BC043504.1</u> , <u>BC008085</u> , <u>BC018494</u> , <u>BX117697</u> , <u>NM_024726.5</u>
UniProt ID:	<u>Q86XH1</u>
Synonyms:	4930465P12Rik; DRC11; IQCA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)