

Product datasheet for **HP215007**

EHMT1/GLP (EHMT1) Human qPCR Primer Pair (NM_024757)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79813
Forward Sequence:	GCTTCAGAAGGTGCTCCTCATG
Reverse Sequence:	CTGAACCAGCATGTGGCAGATG
ACCN:	<u>NM_024757</u> , <u>NM_024757.1</u> , <u>NM_024757.2</u> , <u>NM_024757.3</u> , <u>NM_024757.4</u> , <u>BC047504</u> , <u>BC011608</u> , <u>BC025772</u> , <u>BM681619</u> , <u>BQ028833</u> , <u>NM_024757.5</u>
UniProt ID:	<u>Q9H9B1</u>
Synonyms:	EHMT1-IT1; Eu-HMTase1; EUHMTASE1; FP13812; GLP; GLP1; KLEFS1; KMT1D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)