

Product datasheet for **HP215011**

MOBKL2B (MOB3B) Human qPCR Primer Pair (NM_024761)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79817
Forward Sequence:	TGAGGTTCAAGATCAACAACGAGG
Reverse Sequence:	GGTCGAAGTGGTGGATATAGACG
ACCN:	<u>NM_024761</u> , <u>NM_024761.1</u> , <u>NM_024761.2</u> , <u>NM_024761.3</u> , <u>NM_024761.4</u> , <u>BC033027</u> , <u>BC033027.1</u> , <u>BF511468</u> , <u>BI860336</u> , <u>NM_024761.5</u>
UniProt ID:	<u>Q86TA1</u>
Synonyms:	C9orf35; MOB1D; MOBKL2B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)