

Product datasheet for **HP215015**

CLNMT (CAMKMT) Human qPCR Primer Pair (NM_024766)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79823
Forward Sequence:	CTGAAGAGGAGGTTGGTGCATG
Reverse Sequence:	CCTGTATTGTCAAAGCTGGTAAGG
ACCN:	BC029359 , NM_024766 , NM_024766.1 , NM_024766.2 , NM_024766.3 , NM_024766.4 , BC053733 , BC053733.1 , BC022192 , NM_024766.5
UniProt ID:	Q7Z624
Synonyms:	C2orf34; Cam; CaM KMT; CLNMT; KMT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)