

Product datasheet for **HP215017**

ASAM (CLMP) Human qPCR Primer Pair (NM_024769)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79827
Forward Sequence:	TTACTGGCAGCGAATCCGAGAG
Reverse Sequence:	CTGGTACAGTCCAGAGTAGGAC
ACCN:	NM_024769 , NM_024769.1 , NM_024769.2 , NM_024769.3 , BC009371 , BC009371.2 , BE272517 , BG723094 , BK001245 , NM_024769.5
UniProt ID:	Q9H6B4
Synonyms:	ACAM; ASAM; CSBM; CSBS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)