

Product datasheet for **HP215044**

BORA Human qPCR Primer Pair (NM_024808)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79866
Forward Sequence:	CCTTCGCCTATTTCTTCACCCAC
Reverse Sequence:	CCAGATGAAGCATCAGGAGAATG
ACCN:	BC025367 , NM_024808 , NM_024808.1 , NM_024808.2 , NM_024808.3 , BC007119 , BC044213 , BC056876 , BM923278 , BX648605 , NM_024808.5
UniProt ID:	Q6PGQ7
Synonyms:	C13orf34
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)