

Product datasheet for **HP215056**

ZC3H14 Human qPCR Primer Pair (NM_024824)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79882
Forward Sequence:	CGACCCAGTAATGGCAGAACTC
Reverse Sequence:	CGACTCTTGTTTGGTGCCACT
ACCN:	BC027607 , NM_024824 , NM_024824.1 , NM_024824.2 , NM_024824.3 , NM_024824.4 , BC011793 , BC011793.1 , BC023641 , BC034274 , BX248265 , BX248281 , BX248769 , BX640716 , BX640876 , NM_024824.5
UniProt ID:	Q6PJT7
Synonyms:	MRT56; MSUT-2; NY-REN-37; SUT2; UKp68
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)