

Product datasheet for **HP215067**

ZNF672 Human qPCR Primer Pair (NM_024836)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79894
Forward Sequence:	TTCCGTTGCCTAGAATGCGGTG
Reverse Sequence:	TGTGCGCCCAAGTGTAGGTCAA
ACCN:	NM_024836 , NM_024836.1 , NM_024836.2 , BC068506 , BC068506.1 , BC007385 , BC008872 , BC029310 , BC035140 , BC050348 , BX449043 , NM_024836.3
UniProt ID:	Q499Z4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)