

## Product datasheet for **HP215147**

### PHC3 Human qPCR Primer Pair (NM\_024947)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80012
Forward Sequence:	CCTACAAGTGCAACCACCAGCA
Reverse Sequence:	TGCTGGAGACAAAGTGGGTGGT
ACCN:	<a href="#">NM_024947</a> , <a href="#">NM_024947.1</a> , <a href="#">NM_024947.2</a> , <a href="#">NM_024947.3</a> , <a href="#">BC131772</a> , <a href="#">BC022325</a> , <a href="#">BC026001</a> , <a href="#">BC070164</a> , <a href="#">BC112944</a> , <a href="#">BC131773</a> , <a href="#">BP200133</a> , <a href="#">BX647868</a> , <a href="#">NM_024947.4</a>
UniProt ID:	<a href="#">Q8NDX5</a>
Synonyms:	EDR3; HPH3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)