

Product datasheet for **HP215188**

MZT2B Human qPCR Primer Pair (NM_025029)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80097
Forward Sequence:	AGGTGCTGAGCACCGAGGAGA
Reverse Sequence:	GCTTCAGCAGGTCCACCAGGA
ACCN:	BC017694 , NM_025029 , NM_025029.1 , NM_025029.2 , NM_025029.3 , NM_025029.4 , BC066296 , BC066296.1 , BC014255 , BG036944 , BI546951 , NM_025029.5
UniProt ID:	Q6NZ67
Synonyms:	FAM128B; MOZART2B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)