

Product datasheet for **HP215206**

LRR8E Human qPCR Primer Pair (NM_025061)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80131
Forward Sequence:	GCCGAGTACCTCACCGTGGC
Reverse Sequence:	GCCTCTGATAAGTTCTCCTGGAG
ACCN:	NM_025061 , NM_025061.2 , NM_025061.3 , NM_025061.4 , NM_025061.5 , BC070089 , BC070089.1 , BC022216 , BC108252 , BX538180 , NM_025061.6
UniProt ID:	Q6NSJ5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)