

Product datasheet for **HP215263**

Nucleotide binding protein like (NUBPL) Human qPCR Primer Pair (NM_025152)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80224
Forward Sequence:	CACCAGGAAGTGGAGATGTGCA
Reverse Sequence:	CTCTGCGAAACATCTCAGCACC
ACCN:	NM_025152 , NM_025152.1 , NM_025152.2 , BC024919 , BM978769 , BP274976 , BQ649741 , BX248028 , BX648419
UniProt ID:	Q8TB37
Synonyms:	C14orf127; huInd1; IND1; MC1DN21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)