

## Product datasheet for **HP215305**

### PNPLA3 Human qPCR Primer Pair (NM\_025225)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80339
Forward Sequence:	CCTGTGGAATCTGCCATTGCGA
Reverse Sequence:	GAGCAGACACATCAGCACTCGA
ACCN:	<u>NM_025225, NM_025225.1, NM_025225.2, BC065195, BC065195.1, BC014449, NM_025225.3</u>
UniProt ID:	<u>Q9NST1</u>
Synonyms:	ADPN; C22orf20; iPLA(2)epsilon
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)