

Product datasheet for **HP215307**

TRAF3IP3 Human qPCR Primer Pair (NM_025228)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80342
Forward Sequence:	GCCAGTATGAGGCTCTGAAGGA
Reverse Sequence:	GCCTGGTTCATCTGTAAGTCCC
ACCN:	NM_025228 , NM_025228.1 , NM_025228.2 , NM_025228.3 , BC110302 , BC130417 , BC144139 , BQ277340 , BX648450
UniProt ID:	Q9Y228
Synonyms:	T3JAM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)