

Product datasheet for **HP215321**

ACAD10 Human qPCR Primer Pair (NM_025247)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80724
Forward Sequence:	GCCTCTTCAGATGCCACCAACA
Reverse Sequence:	CACAGAGTTGGCAACGAGGATC
ACCN:	<u>BC015056</u> , <u>NM_025247</u> , <u>NM_025247.1</u> , <u>NM_025247.2</u> , <u>NM_025247.3</u> , <u>NM_025247.4</u> , <u>NM_025247.5</u> , <u>BC015056.1</u> , <u>BC126358</u> , <u>BC003698</u> , <u>BC054886</u> , <u>NM_025247.6</u>
UniProt ID:	<u>Q6JQN1</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)