

## Product datasheet for **HP215322**

### SNAP25 Interacting Protein (SRCIN1) Human qPCR Primer Pair (NM\_025248)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 80725  |
| Forward Sequence: | TGCGCTATCTCAACGACGAGGA   |
| Reverse Sequence: | TTCAGCACCAGTGCCTTCTCT  |
| ACCN:             | <a href="#">NM_025248</a> , <a href="#">NM_025248.1</a> , <a href="#">NM_025248.2</a> , <a href="#">BC029278</a> , <a href="#">BC033233</a> , <a href="#">BC156056</a> , <a href="#">BC172455</a> , <a href="#">BX648057</a>   |
| UniProt ID:       | <a href="#">Q9C0H9</a>   |
| Synonyms:         | P140; SNIP   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)