

Product datasheet for **HP215323**

IQCN Human qPCR Primer Pair (NM_025249)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80726
Forward Sequence:	GGAGACAACGGAGCCACACGT
Reverse Sequence:	CATGCGCGTTCCACGATTGG
ACCN:	<u>NM_025249</u> , <u>NM_025249.1</u> , <u>NM_025249.2</u> , <u>NM_025249.3</u> , <u>BC112932</u> , <u>BC034327</u> , <u>BC127099</u> , <u>BC127100</u> , <u>BQ006839</u> , <u>NM_025249.4</u>
UniProt ID:	<u>Q9H0B3</u>
Synonyms:	KIAA1683
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)