

Product datasheet for **HP215332**

TSEN2 Human qPCR Primer Pair (NM_025265)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80746
Forward Sequence:	GCCTTTTCTTGGTCTATGCTCTG
Reverse Sequence:	GCCATGTAGGTGGTTCTGAACG
ACCN:	NM_025265 , NM_025265.1 , NM_025265.2 , NM_025265.3 , BC019582 , BC019582.1 , BC004178 , BC004211 , BC021975 , BE858588 , BM671402 , BX641631 , NM_025265.4
UniProt ID:	Q8NCE0
Synonyms:	PCH2B; SEN2; SEN2L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)