

Product datasheet for **HP215349**

CYP2S1 Human qPCR Primer Pair (NM_030622)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29785
Forward Sequence:	GATGGACGGTTCAGGAAGCATG
Reverse Sequence:	GGAGAAGGCTTGTAGGATGGTG
ACCN:	NM_030622 , NM_030622.1 , NM_030622.2 , NM_030622.3 , NM_030622.4 , NM_030622.5 , NM_030622.6 , NM_030622.7 , BC033691 , BC033691.1 , BP284240 , BX422464 , NM_030622.8
UniProt ID:	Q96SQ9
Synonyms:	CYPIIS1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)