

## Product datasheet for **HP215354**

### CPEB4 Human qPCR Primer Pair (NM\_030627)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80315
Forward Sequence:	GGATGGTTCACAGCCACTTGAC
Reverse Sequence:	CACACCTCCATATAGCCGATCC
ACCN:	<a href="#">NM_030627</a> , <a href="#">NM_030627.1</a> , <a href="#">NM_030627.2</a> , <a href="#">NM_030627.3</a> , <a href="#">BC036899</a> , <a href="#">BC117150</a> , <a href="#">BC143958</a> , <a href="#">BC143959</a> , <a href="#">BQ061431</a> , <a href="#">BU629646</a> , <a href="#">BX401008</a> , <a href="#">BX538213</a> , <a href="#">NM_030627.4</a>
UniProt ID:	<a href="#">Q17RY0</a>
Synonyms:	CPE-BP4; hCPEB-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)