

Product datasheet for **HP215379**

SEC13L1 (SEC13) Human qPCR Primer Pair (NM_030673)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6396
Forward Sequence:	ACAACGCTCACACCATTTGGCTG
Reverse Sequence:	GTTGTCACAGCCACCTGATGCA
ACCN:	<u>NM_030673</u> , <u>NM_030673.2</u> , <u>BC091506</u> , <u>BC002634</u> , <u>BC006167</u> , <u>BE621430</u>
UniProt ID:	<u>P55735</u>
Synonyms:	D3S1231E; npp-20; SEC13L1; SEC13R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)