

Product datasheet for **HP215395**

ILKAP Human qPCR Primer Pair (NM_030768)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80895
Forward Sequence:	CCGTGAAGAGATGCCTTTGGAC
Reverse Sequence:	GTTGTCTACAGCCAGAACACACG
ACCN:	<u>NM_030768</u> , <u>NM_030768.1</u> , <u>NM_030768.2</u> , <u>BC006576</u> , <u>BC035746</u> , <u>BC050422</u> , <u>NM_030768.3</u>
UniProt ID:	<u>Q9H0C8</u>
Synonyms:	ILKAP2; ILKAP3; PP2C-DELTA; PP2CD; PPM1O
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)