

Product datasheet for **HP215403**

Z DNA binding protein (ZBP1) Human qPCR Primer Pair (NM_030776)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81030
Forward Sequence:	GCAAACCTCCGAAGCCATCCAGA
Reverse Sequence:	CCAAGTTGAGGAATCACCTGGTG
ACCN:	BC028218 , NM_030776 , NM_030776.1 , NM_030776.2 , BC028218.1 , BC131706 , BF511589 , NM_030776.3
UniProt ID:	Q9H171
Synonyms:	C20orf183; DAI; DLM-1; DLM1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)