

Product datasheet for **HP215405**

KCNH6 Human qPCR Primer Pair (NM_030779)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 81033 |
| Forward Sequence: | CGAGAAGGTCTTCTCCATCTGC |
| Reverse Sequence: | ACGCAGCATCTGCGTGTGGTAG |
| ACCN: | BC006334 , NM_030779 , NM_030779.1 , NM_030779.2 , NM_030779.3 , BC006334.2 , BE465143 , BM561167 , NM_030779.4 |
| UniProt ID: | Q9H252 |
| Synonyms: | ERG-2; ERG2; hERG-2; HERG2; Kv11.2 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)