

## Product datasheet for **HP215425**

### **C15orf44 (VWA9) Human qPCR Primer Pair (NM\_030800)**

#### **Product data:**

Product Type:	qPCR Primer Pairs
Gene ID:	81556
Forward Sequence:	TTGGTAGAGGGTCACTGCGACA
Reverse Sequence:	GGCATTCCAAGGAATCGGTGCT
ACCN:	<u><a href="#">NM_030800</a></u> , <u><a href="#">NM_030800.1</a></u> , <u><a href="#">NM_030800.2</a></u> , <u><a href="#">BC007991</a></u> , <u><a href="#">BC007991.2</a></u> , <u><a href="#">BC037823</a></u> , <u><a href="#">BM558754</a></u> , <u><a href="#">BP315446</a></u>
UniProt ID:	<u><a href="#">Q96SY0</a></u>
Synonyms:	DKFZp564O1664
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)