

Product datasheet for **HP215431**

MRPS26 Human qPCR Primer Pair (NM_030811)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64949
Forward Sequence:	TGCGGAGTTCTTCGTGCTGATG
Reverse Sequence:	CACCTTCCTCTGCACCTCGGA
ACCN:	<u>NM_030811</u> , <u>NM_030811.1</u> , <u>NM_030811.2</u> , <u>NM_030811.3</u> , <u>BC013018</u> , <u>BC013018.2</u> , <u>BI115733</u> , <u>NM_030811.4</u>
UniProt ID:	<u>Q9BYN8</u>
Synonyms:	C20orf193; dj534B8.3; GI008; MRP-S13; MRP-S26; MRPS13; NY-BR-87; RPMS13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)