

Product datasheet for **HP215439**

GFOD2 Human qPCR Primer Pair (NM_030819)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81577
Forward Sequence:	GTGGATCTGGTGCATCAGCA
Reverse Sequence:	CATCCACCGATGTTGCTGCCTT
ACCN:	NM_030819 , NM_030819.1 , NM_030819.2 , NM_030819.3 , BC107585 , BC107585.1 , BC000757 , BC002943 , BI826603 , BQ003788 , BQ437486 , NM_030819.4
UniProt ID:	Q3B7J2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)