

Product datasheet for **HP215494**

ADAMTS10 Human qPCR Primer Pair (NM_030957)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81794
Forward Sequence:	TCGTAACCTCGCCTCATCCTGCT
Reverse Sequence:	GTGGTTCACGATGGATTCTGCC
ACCN:	NM_030957 , NM_030957.1 , NM_030957.2 , NM_030957.3 , BC029308 , BC169218 , BF510205 , BF588726 , BI822988 , BM928283 , BX461830 , NM_030957.4
UniProt ID:	Q9H324
Synonyms:	ADAM-TS10; ADAMTS-10; WMS; WMS1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)