

Product datasheet for **HP215495**

SLCO5A1 Human qPCR Primer Pair (NM_030958)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81796
Forward Sequence:	CGAGTATGAGCCAGTCTGTGGA
Reverse Sequence:	CTTGCGACTTTGGACACAGGT
ACCN:	<u>NM_030958</u> , <u>NM_030958.1</u> , <u>NM_030958.2</u> , <u>BC139755</u> , <u>BC016067</u> , <u>BC137424</u> , <u>BC144618</u> , <u>NM_030958.3</u>
UniProt ID:	<u>Q9H2Y9</u>
Synonyms:	OATP-J; OATP-RP4; OATP5A1; OATPJ; OATPRP4; SLC21A15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)