

## Product datasheet for **HP215548**

### SLC25A31 Human qPCR Primer Pair (NM\_031291)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83447
Forward Sequence:	CGGTACAAAGGCATGGTGGACT
Reverse Sequence:	CTCCAGACATGAATAGCTGCTTG
ACCN:	<u><a href="#">NM_031291</a></u> , <u><a href="#">NM_031291.1</a></u> , <u><a href="#">NM_031291.2</a></u> , <u><a href="#">NM_031291.3</a></u> , <u><a href="#">BC022032</a></u> , <u><a href="#">BG192958</a></u> , <u><a href="#">BG773445</a></u> , <u><a href="#">BM563946</a></u> , <u><a href="#">BQ438347</a></u> , <u><a href="#">BU853653</a></u> , <u><a href="#">NM_031291.4</a></u>
UniProt ID:	<u><a href="#">Q9H0C2</a></u>
Synonyms:	AAC4; ANT 4; ANT4; SFEC35kDa
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)