

Product datasheet for **HP215555**

MAD3 (MXD3) Human qPCR Primer Pair (NM_031300)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83463
Forward Sequence:	GAGCATGGTTATGCGTCCCTGT
Reverse Sequence:	CCTGCGCTTCTCCAGTTCATTG
ACCN:	NM_031300 , NM_031300.1 , NM_031300.2 , NM_031300.3 , BC000745 , BC000745.2 , BC032586 , BC041690 , BM849150 , BU620484 , NM_031300.4
UniProt ID:	Q9BW11
Synonyms:	BHLHC13; MAD3; MYX
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)