

Product datasheet for **HP215560**

PUS3 Human qPCR Primer Pair (NM_031307)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83480
Forward Sequence:	CTATCACCGATGTGGGAGAACAG
Reverse Sequence:	TAAAGTCCTCGGAATCCCTGCC
ACCN:	<u>NM_031307</u> , <u>NM_031307.1</u> , <u>NM_031307.2</u> , <u>NM_031307.3</u> , <u>BC004822</u> , <u>BC004822.1</u> , <u>BC013427</u> , <u>BM786445</u>
UniProt ID:	<u>Q9BZE2</u>
Synonyms:	2610020J05Rik; FKSG32; MRT55; NEDMIGS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)