

Product datasheet for **HP215565**

ARFBP1 (HUWE1) Human qPCR Primer Pair (NM_031407)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10075
Forward Sequence:	CGGCTTTTCCTGAAGAAGGGAC
Reverse Sequence:	GTTTCAAAGGCTTCAGAGCAGC
ACCN:	NM_031407 , NM_031407.1 , NM_031407.2 , NM_031407.3 , NM_031407.4 , NM_031407.6 , BC002602 , BC018750 , BC038184 , BC063505 , BC072421 , BC107576 , BG334536 , BU619074
UniProt ID:	Q7Z6Z7
Synonyms:	ARF-BP1; HECTH9; HSPC272; Ib772; LASU1; MRXST; MULE; URE-B1; UREB1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)