

Product datasheet for **HP215587**

SOX7 Human qPCR Primer Pair (NM_031439)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83595
Forward Sequence:	TGAACGCCTTCATGGTTGGGC
Reverse Sequence:	GTCAGCGCTTCCACGACTTTC
ACCN:	NM_031439 , NM_031439.1 , NM_031439.2 , NM_031439.3 , BC071947 , BC071947.1 , BC004299 , BE671627 , BG740632 , BT006693 , BU187815 , NM_031439.4
UniProt ID:	Q9BT81
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)