

## Product datasheet for **HP215594**

### ZMIZ2 Human qPCR Primer Pair (NM\_031449)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83637
Forward Sequence:	CCTGGCTGTAAGCAACCATGTC
Reverse Sequence:	GCCAGTTGGTGTTTCATCTGCCG
ACCN:	<u><a href="#">NM_031449</a></u> , <u><a href="#">NM_031449.1</a></u> , <u><a href="#">NM_031449.2</a></u> , <u><a href="#">NM_031449.3</a></u> , <u><a href="#">BC008049</a></u> , <u><a href="#">BC016968</a></u> , <u><a href="#">BC021924</a></u> , <u><a href="#">BC110435</a></u> , <u><a href="#">BC131780</a></u> , <u><a href="#">BI115680</a></u> , <u><a href="#">NM_031449.4</a></u>
UniProt ID:	<u><a href="#">Q8NF64</a></u>
Synonyms:	hZIMP7; NET27; TRAFIP20; ZIMP7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)