

Product datasheet for **HP215613**

SH3BGRL2 Human qPCR Primer Pair (NM_031469)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83699
Forward Sequence:	GAAGAACAGAGGCAATGGATGTAC
Reverse Sequence:	GACTGTGTTGCTTTCCTTGATTG
ACCN:	NM_031469 , NM_031469.1 , NM_031469.2 , NM_031469.3 , BC040489 , BC040489.1 , BC052987 , BC060799 , BC070059 , BC109043 , BC109044 , NM_031469.4
UniProt ID:	Q9UJC5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)