

Product datasheet for **HP215617**

Espin (ESPN) Human qPCR Primer Pair (NM_031475)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83715
Forward Sequence:	TGCAGACCAAGAACAACCTCCGC
Reverse Sequence:	TTGCGGCTGGAGTCCTGCCTC
ACCN:	<u>NM_031475</u> , <u>NM_031475.1</u> , <u>NM_031475.2</u> , <u>BC156375</u> , <u>BC172457</u> , <u>BM713713</u> , <u>NM_031475.3</u>
UniProt ID:	<u>B1AK53</u>
Synonyms:	DFNB36; LP2654; USH1M
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)