

Product datasheet for **HP215622**

RIOK1 Human qPCR Primer Pair (NM_031480)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83732
Forward Sequence:	CACAGCAGAGATACCATGTCCAG
Reverse Sequence:	ACAACTCCCGAGCCTTGGATTC
ACCN:	<u>NM_031480</u> , <u>NM_031480.1</u> , <u>BC006104</u> , <u>BC006104.1</u> , <u>BC000811</u>
UniProt ID:	<u>Q9BRS2</u>
Synonyms:	AD034; bA288G3.1; RIO1; RRP10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)