

Product datasheet for **HP215641**

FTHL17 Human qPCR Primer Pair (NM_031894)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	53940
Forward Sequence:	ACTTCCTGGAGAGCCACTACCT
Reverse Sequence:	GGTGAGCTTGTCGAACAGGTAC
ACCN:	<u>NM_031894</u> , <u>NM_031894.1</u> , <u>BC069069</u> , <u>BC069538</u> , <u>BC100768</u> , <u>BC100769</u> , <u>BC100770</u> , <u>NM_031894.3</u>
UniProt ID:	<u>Q9BXU8</u>
Synonyms:	CT38
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)