

## Product datasheet for **HP215646**

### GRASP65 (GORASP1) Human qPCR Primer Pair (NM\_031899)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64689
Forward Sequence:	TCATGGAGACTCCTCTTCAGCC
Reverse Sequence:	CACACTGGCATTGCTGTTGTCC
ACCN:	<a href="#">NM_031899</a> , <a href="#">NM_031899.1</a> , <a href="#">NM_031899.2</a> , <a href="#">NM_031899.3</a> , <a href="#">BC075854</a> , <a href="#">BC075854.1</a> , <a href="#">BC008928</a> , <a href="#">BC064624</a> , <a href="#">BC103826</a> , <a href="#">BC103827</a> , <a href="#">BM470731</a>
UniProt ID:	<a href="#">Q9BQQ3</a>
Synonyms:	GOLPH5; GRASP65; P65
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)