

Product datasheet for **HP215650**

FRMD8 Human qPCR Primer Pair (NM_031904)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83786
Forward Sequence:	GTGCACGTCATCGATAGCAGAG
Reverse Sequence:	GAGGAGCTTGTTGACAGGTGTG
ACCN:	<u>BC033851</u> , <u>NM_031904</u> , <u>NM_031904.1</u> , <u>NM_031904.2</u> , <u>NM_031904.3</u> , <u>NM_031904.4</u> , <u>BC033851.1</u> , <u>BC044658</u> , <u>BC044658.1</u> , <u>BC014498</u> , <u>BC051695</u>
UniProt ID:	<u>Q9BZ67</u>
Synonyms:	FKSG44; iTAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)