

Product datasheet for **HP215655**

CTRP6 (C1QTNF6) Human qPCR Primer Pair (NM_031910)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114904
Forward Sequence:	CTGCCTGAGATCAGACCCTACA
Reverse Sequence:	TTGTCACCCTTGCTGCCCTGAG
ACCN:	NM_031910 , NM_031910.2 , NM_031910.3 , BC020551 , BC020551.1 , BM678522 , BM722009 , BQ674594 , BU147265 , BU535692 , BU542015 , BU859813
UniProt ID:	Q9BXI9
Synonyms:	CTFP6; CTRP6; ZACRP6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)