

Product datasheet for **HP215706**

FUCA2 Human qPCR Primer Pair (NM_032020)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2519
Forward Sequence:	ATCAGCCTGAGGTCTGTGGTC
Reverse Sequence:	GCTTACAGATGCTACCAGCTCC
ACCN:	NM_032020 , NM_032020.1 , NM_032020.3 , NM_032020.4 , BC051268 , BC051268.1 , BC003060 , BG719318 , BP214015 , BX537633
UniProt ID:	Q9BTY2
Synonyms:	dj20N2.5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)