

Product datasheet for **HP215713**

TSSK1 (TSSK1B) Human qPCR Primer Pair (NM_032028)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83942
Forward Sequence:	GATGCTGCGTATCCAGAAGGAG
Reverse Sequence:	TGGCTGAGGATCTCGTCGATGT
ACCN:	BC022515 , NM_032028 , NM_032028.1 , NM_032028.2 , NM_032028.3 , BC022515.1
UniProt ID:	Q9BXA7
Synonyms:	FKSG81; SPOGA4; STK22D; TSK1; TSSK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)