

Product datasheet for **HP215718**

Sodium bicarbonate transporter like protein 11 (SLC4A11) Human qPCR Primer Pair (NM_032034)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83959
Forward Sequence:	GCCATTGGTGATTCTGCTGACC
Reverse Sequence:	CGTAGAAGGAGTTGAAGTCCAGG
ACCN:	NM_032034 , NM_032034.1 , NM_032034.2 , NM_032034.3 , BC110540 , BC110541 , BM666654 , NM_032034.4
UniProt ID:	Q8NBS3
Synonyms:	BTR1; CDPD1; CHED; CHED2; dj794I6.2; NABC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)