

Product datasheet for **HP215720**

IFI27L2 Human qPCR Primer Pair (NM_032036)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83982
Forward Sequence:	ATGATGTCCGCAGCAGCCATTG
Reverse Sequence:	CCAACACTGACCCAACAGAGGC
ACCN:	NM_032036 , NM_032036.1 , NM_032036.2 , BC032626 , BC032626.1 , BC022800 , BE737110 , BN000226 , NM_032036.3
UniProt ID:	Q9H2X8
Synonyms:	FAM14A; ISG12B; TLH29
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)