

## Product datasheet for **HP215725**

### FAM172A Human qPCR Primer Pair (NM\_032042)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 83989  |
| Forward Sequence: | GGCATCAAGAAGCTGGCAAAACG  |
| Reverse Sequence: | GGCAATCAGGTAGCATGGACTC   |
| ACCN:             | <a href="#">NM_032042</a> , <a href="#">NM_032042.1</a> , <a href="#">NM_032042.2</a> , <a href="#">NM_032042.3</a> , <a href="#">NM_032042.4</a> , <a href="#">NM_032042.5</a> , <a href="#">BC020584</a> , <a href="#">BC020584.2</a> , <a href="#">BC012007</a> , <a href="#">BP398089</a> , <a href="#">NM_032042.6</a>                  |
| UniProt ID:       | <a href="#">Q8WUF8</a>   |
| Synonyms:         | C5orf21; Toupee  |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)