

Product datasheet for **HP215736**

GAJ (MND1) Human qPCR Primer Pair (NM_032117)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84057
Forward Sequence:	CTTCGAGACCAAAGGGAACAGC
Reverse Sequence:	TCAGTCCATCTGTTAGCAGCTTC
ACCN:	NM_032117 , NM_032117.1 , NM_032117.2 , NM_032117.3 , BC032142 , BC032142.2 , BF975093 , BG613292 , BG758047 , BM820228 , BM833450 , BU533354 , NM_032117.4
UniProt ID:	Q9BWT6
Synonyms:	GAJ
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)