

Product datasheet for **HP215738**

MASS1 (ADGRV1) Human qPCR Primer Pair (NM_032119)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84059
Forward Sequence:	CGAAGAAGCCAATAGGACCGTG
Reverse Sequence:	CCTCATGCCAAAGGCTGTTTCAG
ACCN:	<u>NM_032119</u> , <u>NM_032119.1</u> , <u>NM_032119.2</u> , <u>NM_032119.3</u> , <u>BC034748</u> , <u>BI758469</u> , <u>BX493741</u>
UniProt ID:	<u>Q8WXG9</u>
Synonyms:	FEB4; GPR98; MASS1; USH2B; USH2C; VLGR1; VLGR1b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)