

Product datasheet for **HP215772**

WDR75 Human qPCR Primer Pair (NM_032168)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84128
Forward Sequence:	GCGGCAAGAAAAGGAACTGAGC
Reverse Sequence:	GCTGTAATGCAGTCTTCGTGTGG
ACCN:	<u>NM_032168</u> , <u>NM_032168.1</u> , <u>NM_032168.2</u> , <u>BC040567</u> , <u>BC040567.1</u> , <u>BC006816</u> , <u>BC018689</u> , <u>BI460450</u> , <u>BU619995</u>
UniProt ID:	<u>Q8IWA0</u>
Synonyms:	NET16; UTP17
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)